



# Hyperion™

## ASDO - With Hole

with Hub



### Customised Product\*

#### Features

- 2 asymmetrical occluder discs with hole
- 72 pre-oxidised Nitinol wires
- In-disc and in-waist PET patch with hole
- With hub design
- Screw connection

#### Sizing

- 18 to 42 mm (with Hub)

#### Sizing Hole

- 4 to 10 mm (size dependable)

#### Benefits

- Pulmonary Pressure reduction
- Functional elasticity
- Enhanced shape memory
- Easy positioning
- Anatomical deployment
- Safe attachment/detachment
- High bio-compatibility
- Fast endothelialisation
- Minimised nickel release

#### Delivery

The Hyperion™ ASD Occluder is compatible with the Orbit-A Delivery System

(A) : Left Atrium (LA) disc Ø

(B) : Right Atrium (RA) disc Ø

(F) : Fenestrated hole

- Ⓐ : Left Atrium (LA) disc Ø
- Ⓑ : Right Atrium (RA) disc Ø
- Ⓕ : Fenestrated hole

