

# Hyperion™

## Snare



### Product Information

#### Features

- Nitinol shaft
- Nitinol / Gold Tungsten coil loop wire
- Shape memory

#### Sizing

- L: 129 cm
- Ø Shaft: 0.8 mm / 0.031"
- Ø Snare circle: 15 and 20 mm
- Angle: 90°


#### Benefits

Easy retrieval of atraumatic foreign bodies in coronary and peripheral cardiovascular system and the extra-cranial neurovascular anatomy


- Functional elasticity
- Anti kinking
- High radiopacity
- ≥ 0.040" compatible

#### Delivery

Compatible with catheter systems with a minimal ID of 0.040"



Hyperion™ Snare catching a 0.035" guide wire



Ⓛ : Length of shaft

ⓓ : Ø of snare circle

Ⓣ : Ø of snare shaft

Ⓐ : 90°

# Ordering Information

## Hyperion™ Snare

| Catalogue #            | C-SN-15      | C-SN-20      |
|------------------------|--------------|--------------|
| Effective Length (cm)  | 129          | 129          |
| Ø of Snare Circle (mm) | 15           | 20           |
| Angle °                | 90           | 90           |
| Ø Shaft (mm/inch)      | 0.8 / 0.031" | 0.8 / 0.031" |

- Ⓛ : Length of shaft
- ⓓ : Ø of snare circle
- Ⓣ : Ø of snare shaft
- Ⓐ : 90°

